

GK5-glycerol kinase 5 (putative) Gene

Reference: PTXBC032470

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP:glycerol 3-phosphotransferase 5; GK 5; glycerokinase 5; glycerol kinase 5 (putative)		

For research use only.

GENE DESCRIPTION	glycerol kinase 5 (putative)
GENE ID	256356
ORIGIN SPECIES	Human
ACCESSIONS	BC032470
VECTOR	pDONR223
NCBI SYMBOL	GK5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggggctgctcacggaccggagcagagagcgcaggagccgggtacccggctcgtgctggggctggatgtggcagttctgtgatccgctgccacgtctatgaccggggcgccgggtctgctgctccagcgtgcagaaggtagaaaatttatctcaaatggctggg
tagaaaattgatcctgatgttcttggattcaattgttgcgtataaaagaagcagtcgaagctgcaggatacagatgaatcaattgttggcttggcattcaacacagagagcaactttattacgtggacaagaaaacaggaatcatttcacaactttataagttggcaagacttaagag
ctgttgaaattgtaaaattctggaaattcttcttatgaagataatttcacagttcttccgagtgcttcttcttctagaaagtaaacgactttttacagccagtttggcttcttcaaccacagcagacttcttgagattggctctggattttacagaacttgactgaggtgcaaaaggcagttgaa
gaagaaaattgtcgttgggactattgatacctggtgttatataagctcacaagaggtctgtatatgccacagattttcaaatgctagtacaactggacttttgacccatataagatgtgttggagtgaggatgattacctctctaatttcataccactttctcctacctcctgtgagggacaca
agccacaattttgatcagtgatgaagagatattgggtgctctataccaatagttgcctgttaaggtccctggatatgaccagaatatctgctatatattttgggaaggggaccattgagccagttattctatcatggaagtactttataa