



Reference: PTXBC032485

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIFL; apoptosis-inducing factor 3; apoptosis-inducing factor like; apoptosis-inducing factor, mitochondrion-associated, 3; apoptosis inducing factor, mitochondria associated 3		

For research use only.

GENE DESCRIPTION	apoptosis-inducing factor, mitochondrion-associated, 3
GENE ID	150209
ORIGIN SPECIES	Human
ACCESSIONS	BC032485
VECTOR	pDONR223
NCBI SYMBOL	AIFM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]