



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HXA1-homeobox A1 Gene

Reference: PTXBC032547

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BSAS; HOX1; HOX1F; homeobox protein Hox-A1; HOX A1 homeodomain protein; Hox 1.6-like protein; homeo box A1; homeobox 1F; homeobox protein Hox-1F; lab-like protein; homeobox A1		

For research use only.

GENE DESCRIPTION	homeobox A1
GENE ID	3198
ORIGIN SPECIES	Human
ACCESSIONS	BC032547
VECTOR	pDONR223
NCBI SYMBOL	HOXA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaatgcagaagaaagaactctcttcggaattaccatactctagcagtgaggcagctcgggacgtctcagcccgagctacccttcggaccataggattacaactcttcagctcgtgcgcggtcagcgcccaacagttcggcggcgacgaccgtcttactagtgggcagggggggca
gatcgtcttcgccaccaccaccaccaccaccaccaccaccaccggccgcgtacctactacagactctcgggaactcgtcttcgtctctctccactcaagttgttcggcgaacatgatccagaaacttcagtcagaaacttcagtcgcgtccctacgcccctacggtcttaaacacaggaagacagcagaagt
gtgggttaccctcccgatcgcgtcccgctctgtttactctggaatctctctatctctccatggctcagcatcaccaccaccacagggtttatcgtggggcgcggttgcctcaataactgacctacatgatgacagcgagcacagcgtccgctccgtccgtactataaataactctctgc
ctctctccacgccagccacaagaagcctctgctcgtccccgcctcggagacactcttccagcgcagacttttgactggatgaaggtcaaaagaaacccctccaaaacagggaagtgtggagatgactcagctctgggtcaccccaacgcgtgcgcgaacactctactaccaagcagct
cagggaataacaggaaggttccactctcaacagactcagcagcgccgcgcgcagctgtggagatcgtctcgtatcctccactcagcagagacccaaggtgaagatcgtgttccagagaccgccgaatgaagcaaaaagaacgtgagaagaagggtctgtcccatctctccgccacccc
ccgcaagaacacagagaagccgcgaagaaactctggaactgaagctcctccctcgtctctcccccggctctcactcagacactctgactactccactga

More info: Please contact us for any question: contact@proteogenix.fr