

CLASP1-cytoplasmic linker associated protein 1 Gene

Reference: PTXBC032563

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAST1; CLIP-associating protein 1; multiple asters 1; multiple asters homolog 1; protein Orbit homolog 1; cytoplasmic linker associated protein 1		

For research use only.

GENE DESCRIPTION	cytoplasmic linker associated protein 1
GENE ID	23332
ORIGIN SPECIES	Human
ACCESSIONS	BC032563
VECTOR	pDONR223
NCBI SYMBOL	CLASP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagactctggaagccacaaagactccataagggtggtgagagcggctgaggaggctgctccacactggccagttccatccacccggagcagtgcatcaaggctgctgccccatccagacggcggactacccatcaaccttgctgccatcaagatgcagaccaaag
tcgtcgagaggatcgaaaggagtcattgctgcagctcctgtcgacatcatccaggcttgcctgcagggttatgacaacaccgaaagtagtgcgttaaggccagcgtgttttcttagtggcaatttattccgtaacggagaagacctgaacctcacctgcacagctcacaggagca
agatgaagctactaaacttatataaagaggccgaccaccaacagcaacagcagctcctcctccgatgtctccagcagcagtaa