

NLGN4Y-neurologin 4, Y-linked Gene

Reference: PTXBC032567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	neurologin-4, Y-linked; neurologin 4, Y-linked		

For research use only.

GENE DESCRIPTION	neurologin 4, Y-linked
GENE ID	22829
ORIGIN SPECIES	Human
ACCESSIONS	BC032567
VECTOR	pDONR223
NCBI SYMBOL	NLGN4Y
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccatctgtttaccaccagtttggatacttgatgacctatgttcaagatcaaaatgaagactgccttacttaaacatctatgtgccatggaagatggaaccaataaagagaatgcagacgatataaccagtaatgacctgtgaagataaagatattcatgaacagaacagt
aagaagcctgttatggctatccatgggggatcttacatggagggaaccggttaacatgatgatggcagcattttgccagctatgggaacgtcatcgttatcaccattaactaccgtctgggaatactaggtatgcaaggagcagcttgggtggagctcaaaaatgtttaattattttaaact
tccttctactaatttaataaatttttttga