

## NHLRC2-NHL repeat containing 2 Gene

Reference: PTXBC032598

|                |                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | NHL repeat-containing protein 2; 1200003G01Rik; novel NHL repeat domain containing protein; NHL repeat containing 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | NHL repeat containing 2                                                                          |
| <b>GENE ID</b>          | 374354                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC032598                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | NHLRC2                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcaggagctcatcagatagggcactcctgctggactctggcaactgccaaagaaatgagtaacaaaaggaaacctgccttaggtttgctggaagtggaaatgaagagaatcgaacaatgcctatcctcacaaggcaggtttgccaaacctcaggccttcttggcctctgaa  
gatccctggagctgctgtttgtagcagatagtgagagcagtagagaaacgtttcactgaaagatggagcagtgagcaccctgtagggagagaagagacccatgaattattttgctttggtgagttgatggagtaggaatcaatgcaaagctcaacacccctggagtaacat  
gggacaaaaaaggaaattactttatgtgtagcagactcctacaatcacaagattaaagtgtggatccaaaaacaaaaactgtacaacattagcaggaactggagacacaaataatgttacagttccagttttacagagtcacatttaatagaaccaggaggtgtgtatggagagaatggag  
aattattatgtagcagacaccaataatcatcaatcaatgaatgagtgatttagaaactaaatggtatctgtgctcccatcttcagatctgaaatgctgtggttagatggccggttccctagtagaaaaacagaagacattacccaaactacataatctgctcgaagcattaggttcccccgtg  
actgctgtgctggccagactcttcagttcaactcagattagacacccatcaggatcaaaagctaaactgaaggatgtagcagttgctgttttctaacagctgaaggcaatgaatggctacttcaaggacagatagcagctggagatagagaacatttcagtcacaaacaatttcaactaca  
aattcctgatgtgcttactcactgaagccattgtatctgtagtggtgtttcttattactgtagtcagacagcagtgctgtgtatgatgaaggcaattttgtcagtcagcctttacaataacgatacacagcaagggtgcatagctccagtagagctcaggtatgtatttag