



Reference: PTXBC032600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALS4; AOA2; SCAR1; bA479K20.2; SEN1 homolog; amyotrophic lateral sclerosis 4 protein; senataxin		

For research use only.

GENE DESCRIPTION	senataxin
GENE ID	23064
ORIGIN SPECIES	Human
ACCESSIONS	BC032600
VECTOR	pDONR223
NCBI SYMBOL	SETX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]