

ATOH7-atonal homolog 7 (Drosophila) Gene

Reference: PTXBC032621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Math5; NCRNA; PHPVAR; RNANC; bHLHa13; protein atonal homolog 7; atonal homolog 7; atonal homolog bHLH transcription factor 7; class A basic helix-loop-helix protein 13; helix-loop-helix protein hATH-5; atonal bHLH transcription factor 7		

For research use only.

GENE DESCRIPTION	atonal homolog 7 (Drosophila)
GENE ID	220202
ORIGIN SPECIES	Human
ACCESSIONS	BC032621
VECTOR	pDONR223
NCBI SYMBOL	ATOH7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtcctgcaagcccagcggcccgccggcgagcgcggttgacccccgtgcggggcgccagcagtgccgggcacgtgcgcggggcgccggcggtggagagcgcgccgagcggcgccgagcggccgcatgaggggctc
aacactgccttcgaccgttacgcagggtgttccccagtggggcccaggataaaaagctgtccaagtacgagaccctgcagatggcctgagctacatcatggctctgacccggatcctggccgagccgagcgaatcggctcggagcgggactgggtgggtctccactgtgagcactt
cgcccgcgaccactacctccgttccggcgcggaagctgcggggcgagagcgagctgtacagccagagactcttcggcttcagcccgagccctccagatggccacctag