



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GNMT-glycine N-methyltransferase Gene

Reference: PTXBC032627

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HEL-S-182mP; glycine N-methyltransferase; epididymis secretory sperm binding protein Li 182mP

For research use only.

GENE DESCRIPTION	glycine N-methyltransferase
GENE ID	27232
ORIGIN SPECIES	Human
ACCESSIONS	BC032627
VECTOR	pDONR223
NCBI SYMBOL	GNMT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggagacgctgtgacggaccgccctctcctgggggtggcgcgcgaaggccctcccgaccagctacgcgcggggagcgcgcgcgctgtggcagctctatctcgagacaccgcagccgcaccgccgagtaacaaggcatggctcgttggctgctgcgcgcagcagcctg
ccagcgggtgctcgcagctagctgttggcactgggggtggactcattatgctgttgggaaggggcttcagtgtgacgagtggttgatgccagtgcacaagatctgaaagtatgcattaaaggcgcgttggaaaccggcgccagcagagcccgcttcgacaagtgggtctatcgaagaagccaac
tggatagctctggacaagaatgctgtgcccgatcagcagaaggggctgttgatgctgtctatcgtccttggaaaagcttgcctacttgcagcagctcaaaagggagcagcagtgagcagcgcctgcctgcgaataacatctcgagatctgcccggcggggcgccctactgtgctatgcatc
gcaactacgcagcaatctcftagtcagcgcgtgcgcacccgggaagaacatctactataaagtgacttgcacaaagcagctacatcagctgctatgctgagcaacaaggcccaatggcgtgacctctgactatagctgtgcaggtggccggggctgcggccaggtgctctctgcgc
ttgagttaagtcctctctactactaccacactctctgcactcttcacggctctcgaagcagctctcgaagtgaaagtcacgacagcgtctctggcgcaactcaagccttcaagcagcgccaaacacatctccctcactattcatccagctctcaagacagacagactca