

C16orf61-chromosome 16 open reading frame 61 Gene

Reference: PTXBC032631

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf61; 2310061C15Rik; DC13; COX assembly mitochondrial protein 2 homolog; C-x(9)-C motif containing 2; C-X9-C motif containing 2		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 61
GENE ID	56942
ORIGIN SPECIES	Human
ACCESSIONS	BC032631
VECTOR	pDONR223
NCBI SYMBOL	C16orf61
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacctgacttattccacactgcacactgaagaatgcaacgtcttgattacactgcttaaggaatgacacaaaacacacattctgaaatttttggtattgtaattgatgtgatcgggagttgagaaatgcctgaagaatgagtagtagaaacaggaccaagagcaggagcatg
gcattgcaatgcgaagaaacttttaacctccagagggaatccgaaaaataa