



Reference: PTXBC032647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPPH; MPPH1; P85B; p85; p85-BETA; phosphatidylinositol 3-kinase regulatory subunit beta; PI3-kinase subunit p85-beta; PI3K regulatory subunit beta; phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; phosphatidylinositol 3-kinase, regulatory subunit		

For research use only.

GENE DESCRIPTION	phosphoinositide-3-kinase, regulatory subunit 2 (beta)
GENE ID	5296
ORIGIN SPECIES	Human
ACCESSIONS	BC032647
VECTOR	pDONR223
NCBI SYMBOL	PIK3R2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr