

PMAIP1-phorbol-12-myristate-13-acetate-induced protein 1 Gene

Reference: PTXBC032663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APR; phorbol-12-myristate-13-acetate-induced protein 1; PMA-induced protein 1; adult T cell leukemia-derived PMA-responsive; immediate-early-response protein APR; protein Noxa		

For research use only.

GENE DESCRIPTION	phorbol-12-myristate-13-acetate-induced protein 1
GENE ID	5366
ORIGIN SPECIES	Human
ACCESSIONS	BC032663
VECTOR	pDONR223
NCBI SYMBOL	PMAIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgggaagaaggcgcgaagaacgtcaaccgagccccgcgggctccagcagcagcgccgggtacggcgggtacggcgaggacgaagcgggattgctgattgggatgcagctgcgtttcaccaggggcaaaaagctcttctctctcttctctctgcccacttgccttccccggggccacgaggaacaaagtgcagtagctgggaagtcgagtgctactcaactcaggagattggagacaaactgaactccggcagaaactctgaatctgatatccaaactctctgctcaggaacctgactgcatacaaaactgcatgaggggactccttcaaaagagttttcaggagggtgcacgtttcatcaattgaagaagactgcattgtaattga