

EFNA1-ephrin-A1 Gene

Reference: PTXBC032698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B61; ECKLG; EFL1; EPLG1; LERK-1; LERK1; TNFAIP4; ephrin-A1; TNF alpha-induced protein 4; eph-related receptor tyrosine kinase ligand 1; immediate early response protein B61; ligand of eph-related kinase 1; tumor necrosis factor, alpha-induced protein 4; e		

For research use only.

GENE DESCRIPTION	ephrin-A1
GENE ID	1942
ORIGIN SPECIES	Human
ACCESSIONS	BC032698
VECTOR	pDONR223
NCBI SYMBOL	EFNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttcctctggcccccctcttgggtctgtgctgcagctgtgccgctgtgatcgccacaccgtcttctggaacagttcaaatcccaagttccggaatgaggactacaccatacatgtgcagctgaatgactacgtggacatcatgtccgcactatgaagatcactctgtggcagacgc
tgccatggagcagtagtactgtacctgtggtgagcatgaggagtagcagctgtgccagccccagtcacaggaccaagtcgctggcagtgcaaccggccccagtgccaagcatggccgggagaagctgtctgagaagttccagcgttcacacctttcaccctgggcaaggagttcaaa
gaaggacacagctactactacatcctcaaacccatccaccagcatgaagaccgctgcttgagggtgaagtgactgtcagtggaataatcactcacagtcctcaggcccatgtcaatccacaggagaagagactgcagcagatgacccagaggtgcgggttctacatagcatgcgtcac
agtctgccccacgctcttcccacttgctgactgtgctcttccacttctgtgctgcaaacccgtga