



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C6orf223-chromosome 6 open reading frame 223 Gene

Reference: PTXBC032706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf223; chromosome 6 open reading frame 223		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 223
GENE ID	221416
ORIGIN SPECIES	Human
ACCESSIONS	BC032706
VECTOR	pDONR223
NCBI SYMBOL	C6orf223
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]