

TRAPPC1-trafficcking protein particle complex 1 Gene

Reference: PTXBC032717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BET5; MUM2; trafficking protein particle complex subunit 1; BET5 homolog; MUM-2; melanoma ubiquitous mutated 2; multiple myeloma protein 2; trafficking protein particle complex 1		

For research use only.

GENE DESCRIPTION	trafficking protein particle complex 1
GENE ID	58485
ORIGIN SPECIES	Human
ACCESSIONS	BC032717
VECTOR	pDONR223
NCBI SYMBOL	TRAPPC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgtccacaacctgtacctgtttgaccggaatggagtggtgtgtgactacagcgaatggcaccgcaagaagcaagcagggattcccaaggaggaggagataagctgatgtacgggatgctctctatccgctcgtttgtcagcaagatgtccccgctagacatgaaggatgcttc
ctggccttccaaactagccgtttacaaactccattactacgagacgcccacfgggatcaaaagtgtgatgaatactgacttggcgtgggacccatccgagatgtgtgcaccacatctacagtgcgtgtatgtggagctgggtggaagaatccccctgtgccgctgggccaactgtgcaa
agtgtgctctttcgtccccgactggactcctatgttgcgtctgtgcccttctctccgcccggcgtggctga