

LOC100125556-family with sequence similarity 86, member A pseudogene Gene

Reference: PTXBC032772

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	family with sequence similarity 86, member A pseudogene
GENE ID	100125556
ORIGIN SPECIES	Human
ACCESSIONS	BC032772
VECTOR	pDONR223
NCBI SYMBOL	LOC100125556
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccggtccttctcagaactcatcaaaaagacgtgctgtattgccagaagccatcggtgtgctggtcgggtcctgctgaggtgctgctgctgccgggagcaccagcgggctcctga