

FABP1-fatty acid binding protein 1, liver Gene

Reference: PTXBC032801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FABPL; L-FABP; fatty acid-binding protein, liver; fatty acid binding protein 1, liver; liver-type fatty acid-binding protein; fatty acid binding protein 1		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 1, liver
GENE ID	2168
ORIGIN SPECIES	Human
ACCESSIONS	BC032801
VECTOR	pDONR223
NCBI SYMBOL	FABP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtttccggcaagtacaaactgcagagccaggaactttgaagccttcataaggaatcggtctgccggaagagctatccagaaggggaaggatcaaggggtgtcggaatcgtgcagaatgggaagcactcaagttaccatcaccgctgggtccaaagtatccaa
aacgaattcacgggtggggaggaatgtgagctggagacaatgacaggggagaaagtcaagacagtggttcagttgggaagtgacaataaactggtgacaactttcaaaacatcaagttctgtgaccgaactcaacggcgacataatcaccaatccatgacattgggtgacattgtctca
agagaatcagcaagagaatttaa