

OCIAD2-OCIA domain containing 2 Gene

Reference: PTXBC032808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OCIA domain-containing protein 2; ovarian carcinoma immunoreactive antigen-like protein; OCIA domain containing 2		

For research use only.

GENE DESCRIPTION	OCIA domain containing 2
GENE ID	132299
ORIGIN SPECIES	Human
ACCESSIONS	BC032808
VECTOR	pDONR223
NCBI SYMBOL	OCIAD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcagcgtctgctgtggaaccaagataaagatgccatgtccaccaccaagcaagcagagcctgtgtttgtccaaaactgcacatccacagagcagagatctcaagattatgcgagaatgtcaggaagaaagttctggaagagagctctgccttttcttctaag
catgctgtcaccaggagactagctaccagggttaattggcagctaattctagatttgatcattgccaaagtgcacgcactgcctccttacctgtgaggaaatgcaaaataa