

RPS29-ribosomal protein S29 Gene

Reference: PTXBC032813

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA13; S29; 40S ribosomal protein S29; ribosomal protein S29		

For research use only.

GENE DESCRIPTION	ribosomal protein S29
GENE ID	6235
ORIGIN SPECIES	Human
ACCESSIONS	BC032813
VECTOR	pDONR223
NCBI SYMBOL	RPS29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcaccagcagctgtactggagccacccgcgaaaatcggccaggggtctctgtgtgtctgttcaaacggcaggtctgatccggaatatggcctcaatgtgccgccagtggttcgtcagtagcgaaggatatcggtttcattaagtggactaa