

KRT222P-keratin 222 pseudogene Gene

Reference: PTXBC032815

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KRT222P; KA21; keratin-like protein KRT222; keratin 222, type II; keratin-222 pseudogene; keratin 222		

For research use only.

GENE DESCRIPTION	keratin 222 pseudogene
GENE ID	125113
ORIGIN SPECIES	Human
ACCESSIONS	BC032815
VECTOR	pDONR223
NCBI SYMBOL	KRT222P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaactgtcccagctactcaatgagatcagggaactatgaaagatcctaccagaaatcagatagagacgggtgctctcaacaaggatccagctagaagaagacataagcaaaaaatggacaagatgaagagctttgaaggcagctcaagcagaactcaaggagcccgac
gccagtggcaccacctgcaagtggaaatgaaatcctccatgctgtggaaggggccttgaactccctacatgccagcggagcagcattaccagatgcagctgcaagacctagagactgtgattgaaggactagaaaaagagctacagggaagtaaggcgcggcatcgaaaaagcagctt
caagagcacgagatcttctcaacacgaagatgaggctggaacaagaatagcaacttatgccacctcctagaaaaagaagaatcagatatattgtgtatccaaaggtgggaaaaagacaaaaagcctaccacaagtagagttgttttttaccatcagccattataatgaaat
ctttactacaaaagtcacaaaagatgagaatgaaatgtagaaacagtaacaaacaggaatccttaaatggaggtatcgttaaggagagcactgaagctcatggcactatcagacagagaaagtggaatgaagttaaaagaatgggaagggtctttctttaagataacccctcattg
aggaaaaagctgtttctctcatttgaatcatttagcagccactgatgaagggtgttagagactaagcaggataatctaccagatatagaagtcaggcttatcatgagaagatcatgcagattccctctatcaaacctccatcaacagctaataa