



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

C20orf144-chromosome 20 open reading frame 144 Gene

Reference: PTXBC032818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf144; dJ63M2.6; bcl-2-like protein from testis; belt; chromosome 20 open reading frame 144		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 144
GENE ID	128864
ORIGIN SPECIES	Human
ACCESSIONS	BC032818
VECTOR	pDONR223
NCBI SYMBOL	C20orf144
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

ORF1 (Only for Reference, product provided without STOP codon):
atggcaaacatctagttcccaacaaaggacacccaagcaccgcccgaaggagagagcgctctgacatggcagcaagcctgttctggaatcctgttctcaaacacctcactcgcgaagaagccggcttaccgagatcgtcgtctatctccccctggacaaagcgcgacccgtctgcgccaa
cgtctgcgcgcgagtgattgacatcgtctcgcgcgtcgtcggctgcggctcccccggccgcacacagcaggtatgcgcggagagcggcggaagctgcgagcgcgagcgcgagatcgcgctgactctcgtctcgtcgtcgcggcgcaagaggcggcgccgccgaagaggcgggcgccgaggcagctctc
gacgtcgcgcgcgcctctctcgtctcgtcgtcgtcgcgcgaagcgtctccccgcgaagccgcgcgcgcgcgcgaagcagcgcgcgcgaacgcacgcgcgcgaacgcctgcgcctccctccccgcgacttttaa