

## TMEM67-transmembrane protein 67 Gene

Reference: PTXBC032835

|                |                                                                                           |                     |                |
|----------------|-------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | JBTS6; MECKELIN; NPHP11; TNEM67; meckel syndrome type 3 protein; transmembrane protein 67 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | transmembrane protein 67                                                                         |
| <b>GENE ID</b>          | 91147                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC032835                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TMEM67                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcgggttggcctctatccgccggcgccgtgaccggttctctgtgttctccctcgtctctacagccagacctctcttcccttccagcagccggagaagtgcgacaacaacagctcttgatatctcgcctctctgtgttctctgtggagctaaccagaggcaagatgc  
ccgaggaaacttcatgtgtatgtctaccaggaatttcagatgatcttaataatggaggacgtctattattgtaaaaagtgccagaaaaatgaaaggtgttacagaagatgctggaactgacttcttgcctagtgacttaactgccgaagaaaaatgtcactgtccattggccatattttagt  
ggaaagagacattaatggaacattgttctcagcaactgtgagctctgctgatggaatgaaactctttatggtagtaaatgctttaggagacaggtgctgctcgtatgtgagccaacttgttaataccagcaggtcctgtgcatgttcagaacccaactttaaagggttatttca  
gcagcacaggggaatttctctacgtagaatttcagctgcacgttatggagaagttggcatgtctttaaactcagaatggttgcgaagtatttgaatcatcagcagctgcatgttgggtatattgccaatctaactctgtcaagctcttggaaatatgtgtgatgaacatgaattctacgacttt  
gccacatttgatgcatgttgactatttcagtttatctttaaataactgctggactgagcactgttcttctatttcttggagacagaacttcttctgctgttattgagaccagtttaggattagcacctcaagtgtcagttctacctctcttctcaaaatttcagttttaaaggagaaaaaccagaa  
tacaaaactgaagtgttggctcttctctatgataagaggaatttctcaagtggcaactttagaaggaggtgttttacagcttggcagacacagagacaaggctaaatgctgcttatttggacaacccaacaaatgtgagattctctatctaaatcttaattgactttccca  
ctctataatttatgatgtgtaccttgaataactgatgaaaaatcaacatcaatataatttggcgtgctgtgtttaaacttaaaacttcaacataataagatatatttgaaccaagacagcagcaactctggaaggtggcttctaactcggcgcaatttcttagtggatgcagtaagtggacgagaaaaatgact  
taggaactcagcaagagtaattcaggttgcactcaaatatcactgagtgccactgttaccacaacataaataatggaacatctaccctccctaatcaccattgectacagtgcattgatatcaaaagatgccaacagccagctgtggaaggttcttctcagtcacatatgaatggatcat  
ggagaagcacatgtccagacagatatgttcttgggtgtatgggtggcgtagctgttttagcatctctttagaacagacagcaggtggaagagcgcatgggagtcctcatgattgattacagacaggtgtgaaattcttgggtactatgtctgtgatctgccaatgtttcttcatcacagtg  
ggaacaggtcttactggcttatttctcaagcacagaagctctgtctgttctgctccaatgccaatcaggaagaacgttctgacttattgttgatgtgcttctgctgaaggcactacaatttgcataagctcatatccagattacaatagatgtattcttattgattgggagcgacacaa  
aggaaaggtcttaagctgttgaaggtgaggtgtgtgtacgaagtgcactgttctgtgaagcatatggagaacataattttagcaaatgaatgaatgaatgaatcagactgtgagaaaaataactcacttcaagtactactgtctctcttcttggaggttgggtgattcaagaacttagca  
ataatggactcatctctagctcttctagaaacccacctagctacatagctctttagctgctttagatagatcagtgctgtctgtcttcttggctagccattggaattatacaggtgctgttcttctgtctttagagagattatagaagataaaatcagacagctgttggattatgtctctatgagta  
accagatgagacaacattatgacagaattcatgacacactaataaggaaaaatggtcctgctagactactgagttcatcagcaagacttctgagcagatataaaagcatatcatatgatgaataaatttctgctccttcatgacatgttcataaggaaatggattactttataaaagataagtt  
gctcttgaagaattcttgaatggaattcatggaaccaatggaaaaaagcatcttttacaatgatgaaggttattcttccagcaggtgctcttattatggaatgaagcactcttctattttagctgtgttcttctgtgttgattggcttcccaatttatttagcactctctacatact  
acaacaagagatttttagatatactcctaatacagtaggacaaaaaattggcatccaaaacttgggtgatcaaaagattttgatttaa

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)