

ISCA2-iron-sulfur cluster assembly 2 homolog (S. cerevisiae) Gene

Reference: PTXBC032893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBLD1; ISA2; MMDS4; c14_5557; iron-sulfur cluster assembly 2 homolog, mitochondrial; HESB-like domain-containing protein 1; iron-sulfur cluster assembly 2		

For research use only.

GENE DESCRIPTION	iron-sulfur cluster assembly 2 homolog (S. cerevisiae)
GENE ID	122961
ORIGIN SPECIES	Human
ACCESSIONS	BC032893
VECTOR	pDONR223
NCBI SYMBOL	ISCA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccgctgggggtcgtccctaacggccgcagcagagagcgggtcactccctggccgaggggagggctcctcacggcctccctgggacccagcgcgtcgggagggcgtcgtcctccagcccgagggcgcgaaggcagatctgcctcacagacagttgcgtccag
aggctttggaaatcaccgaagggtcagaattcctcaggctgcaagtgaggggaggtggtcctccggtccaatacaaatctcactggatcacagttatcaaccccgacgacaggggtatttgaacaggggtgggcaagagtggtggtgactctgatacttgccctcgtgaaaggggc
ccaggtggacttcagccaagaactgatccgaagctcatccaagtgttgaacaatcctcaagcacagcaaggctgctcgtgggtcactcttctctatcaaaattga