



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MOBKL2B-MOB1, Mps One Binder kinase activator-like 2B (yeast) Gene

Reference: PTXBC033027

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBKL2B; C9orf35; MOB1D; MOB kinase activator 3B; MOB1, Mps One Binder kinase activator-like 2B; mob1 homolog 2b; monopolar spindle 1 binding, MOB1, domain containing; mps one binder kinase activator-like 2B		

For research use only.

GENE DESCRIPTION	MOB1, Mps One Binder kinase activator-like 2B (yeast)
GENE ID	79817
ORIGIN SPECIES	Human
ACCESSIONS	BC033027
VECTOR	pDONR223
NCBI SYMBOL	MOBKL2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcacatgccctgaagcaggtattcaacaggcaagacgtctccaccacaagggaaatttgaaactcgcacacagaggtttgagctgcacaacgggctcaggcatccctaaactcgggtgtggacctgaaggcgctgtgcagtgtccacgtggggaggaccagaatgactgggtggcagtcacatgtgggtgaccttcttaactggatcaaccctcatctatggtccacactctgtgagttctgcaccgagcggacctgtcctgtgatgtcagggggcccaaatatgatgtactgtggcgaggtatgctcaagtataaagaccacaacgcctgccagctcccaggtacatgaacctctctctggattgtgaatgtctcatgatcaacaacggaaatattccaaatgcgtggctgtctccctcccaagaactctcttcagatctcagtagaactctgtgcgcctcttctcgggtcttctcggctcagctctatatacccaactctgcagcgggtcattgtgatgtgtcagagggcccatgtcaaca

More info: Please contact us for any question: contact@proteogenix.fr