



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

SERINC1-serine incorporator 1 Gene

Reference: PTXBC033029

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TDE1L; TDE2; TMS-2; TMS2; serine incorporator 1; tumor differentially expressed protein 2		

For research use only.

GENE DESCRIPTION	serine incorporator 1
GENE ID	57515
ORIGIN SPECIES	Human
ACCESSIONS	BC033029
VECTOR	pDONR223
NCBI SYMBOL	SERINC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr