



Reference: PTXBC033030

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	mutS protein homolog 4; hMSH4; mutS homolog 4		
For research use only.			

GENE DESCRIPTION	mutS homolog 4 (E. coli)
GENE ID	4438
ORIGIN SPECIES	Human
ACCESSIONS	BC033030
VECTOR	pDONR223
NCBI SYMBOL	MSH4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr