

C2orf52-chromosome 2 open reading frame 52 Gene

Reference: PTXBC033054

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf52; long intergenic non-protein coding RNA 471		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 52
GENE ID	151477
ORIGIN SPECIES	Human
ACCESSIONS	BC033054
VECTOR	pDONR223
NCBI SYMBOL	C2orf52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagaagctaaggacaatggcagtagagatgaggttctgtccctcataagaattgcaggaagaataccactgtccgggaaagaaggaggagaaagtctctggctcctgtgtttgctgagaagtaataatcaccaagcaggagggagccaagctcaaggaccgtgagagcca
ccaggagaatgaagacaggaacagtgagttggaccagatgaggaagataaagaatcctctgtaggggttcccgatgagtggtgtgagttagagacaagctgtgtgtgtgccattctacagcacttggggagaggttctgttaa