

CIB2-calcium and integrin binding family member 2 Gene

Reference: PTXBC033108

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFNB48; KIP2; USH1J; calcium and integrin-binding family member 2; DNA-dependent protein kinase catalytic subunit-interacting protein 2; calcium and integrin binding family member 2		

For research use only.

GENE DESCRIPTION	calcium and integrin binding family member 2
GENE ID	10518
ORIGIN SPECIES	Human
ACCESSIONS	BC033108
VECTOR	pDONR223
NCBI SYMBOL	CIB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactacaggaagagccccatgctccacgtgcccagcctcatccagatgccagagctccgggagatccctcaagaaagatcgtggcggcggtttccgaggatggtgaggggaacctcacttcaacgactttgtgacatgtttccgtgctctcgagtcggtcccg
agagctcaaggcaactatgcttcaagatctatgacttcaacactgacaacttcatctgcaaggaggacctggagctgacgtgcccggctcactaagtcagagctggatgaggaggaggtggtgcttgtgtgcgacaaggtcattgaggaggctgactggatgggtgacggcaagct
gggctttgctgacttcgaggacatgattgccaaggccctgacttctcagcactttccacatccggaatctga