

ZNF697-zinc finger protein 697 Gene

Reference: PTXBC033126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 697		
For research use only.			

GENE DESCRIPTION	zinc finger protein 697
GENE ID	90874
ORIGIN SPECIES	Human
ACCESSIONS	BC033126
VECTOR	pDONR223
NCBI SYMBOL	ZNF697
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtccgagtcggcgagacctcagcgtcagctgcacctcttcacgcacaagcgcacgcactcgggtgagcggccctacgtgtgccgagtcggggaagggttcgggctgaactcgacactggtgaaccacctgcgcgtgcacaccggcgagaagccctccgtgtggc
cagtcgcgagaagcgttcagcgaactctccacgtcagcagcaccagcgcacgcacacggcgagaagccctacacgtgcatcgagtcgggcaagagctttatccagagctccacctgatccgccaccgcgcacacggcgacaagccgacaagtgtcgggctgcg
gcaaaggctccgctataaaacgcacctcgcgcagcaccagaagctgcacctgtgttag