



Reference: PTXBC033308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein DZIP1L; DZIP2; DAZ interacting protein 1-like; DAZ-interacting protein 1-like protein; DAZ interacting zinc finger protein 1 like		

For research use only.

GENE DESCRIPTION	DAZ interacting protein 1-like
GENE ID	199221
ORIGIN SPECIES	Human
ACCESSIONS	BC033308
VECTOR	pDONR223
NCBI SYMBOL	DZIP1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]