

STX11-syntaxin 11 Gene

Reference: PTXBC033519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FHL4; HLH4; HPLH4; syntaxin-11; syntaxin 11		
For research use only.			

GENE DESCRIPTION	syntaxin 11
GENE ID	8676
ORIGIN SPECIES	Human
ACCESSIONS	BC033519
VECTOR	pDONR223
NCBI SYMBOL	STX11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagaccgctagcagaactctggactgtccaagcaatagaccagcagttccagacggggacgatgagttgactgccccacgaggacatcctggagtcctgtaccgagacatccgggacattcaggatgaaaccagctgctgtggc
cgacgtgaagcggctgggaaagcagaaccccgttctcagctccatgcccgcctcagcagcatcaagcgcgacaccaactccatgccaaaggccatcaagcccgggcgaggtcatccactgcaagctgcgcgccatgaaggagctgagcaggcggctgagggccagc
acggcccgactcggcagtgccgcgcatcttgcgggcgagtaaacgcgtcacctcaccttccagcgcgcatgacgactacaaccagccgagatgaagcagcgcgacaactgcaagatccgcatccagcggcagctggagatcatggcgaaggagctcggcgacca
gatcgaggacatgttcgagcagggtaagtgggacgtgtttccgagaactgctggccgacgtgaaggcgccgcccgcctcaacgagatcgagagccgaccgcgaactgtctgcctggagagccgcatccgcgacgtacacgagctcttctgcagatggcggtgctgtgt
ggagaagcagggcagaccctgaacgtcatcgagctcaacgtacaaaagacggctgactacaccggccaggccaaggcgaggtgcgggaaggccgtgcagtagaggagaagaacccctgccggaccctctgtgcttctgtctccctgctcaagtag