

MMGT1-membrane magnesium transporter 1 Gene

Reference: PTXBC033588

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EMC5; TMEM32; membrane magnesium transporter 1; ER membrane protein complex subunit 5; transmembrane protein 32		

For research use only.

GENE DESCRIPTION	membrane magnesium transporter 1
GENE ID	93380
ORIGIN SPECIES	Human
ACCESSIONS	BC033588
VECTOR	pDONR223
NCBI SYMBOL	MMGT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgtcgtgtggaaggggctggtggcgcgtctctttgcccctagccacgcccctttccgctgcgcagcatcgttcttatatgcgattaacagaaaaagaagatgaatcactgccaatagatatagttctcagacactctggcctttgcagttacgtgtacggtatagttcatat
tgcaggagagtttaagacatggatgccacttcagaaactgaaaaataagacatttgatagcttaaggaaacccatcctttatgtatttaacatcgtgtgcgagttactttccggccttcggatagcgaattctcaaaccaagatgcattgtccttaacacatcattgaagttacgaaaactc
gaatcactgcgtcgttaa