

TEX2-testis expressed 2 Gene

Reference: PTXBC033661

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HT008; TMEM96; testis-expressed sequence 2 protein; testis expressed sequence 2; transmembrane protein 96; testis expressed 2		

For research use only.

GENE DESCRIPTION	testis expressed 2
GENE ID	55852
ORIGIN SPECIES	Human
ACCESSIONS	BC033661
VECTOR	pDONR223
NCBI SYMBOL	TEX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactcagcaaaataaagctcccctactttatgaatgagctcactctgacgggaacttgacatggcgctggtgtgccaataatctccaggccttcaagccttacgttgatcaccaaggactctggattgattggaaatgtcctacaatgggtcctttctgatgactctcgagacaaaatg
aatgtgacaaactaggttaagagcctctgttgaaagcctgaaaggttgagaaattggcaagaaggttgacagcccggtgctctgtgctggacagcgtgaggaatcctccagcgtgctcctccgaggaagacgatgccccagagcccagcggggagacaaacagctcc
tcccagggctgaaggtgacgttgagggtcatgaacaagtaagattatgaggtttgtgataaaattaccaagtcataaatttccaaaaagcaacagagacagagtttattaaaaagaagatcgaagaagctcacaacacccctgctgctcactgttgaaagtacaagaatgtagagaa
ccttggcgggtcaactccaccaccccgactgaccgagtgatggtttccgaaagccaccacatgtggagctgaaagctcggccaaaactggagagagagaagtacttttagttcatgtgacagactggatagagaagaactggagcaagagtttcagaaagttttgtcatgcc
aacatggatgatgtttatatactataatgcactcagccatggacccctgctctacttctgctcctctgaaagaccacctgtggaggtgctgctgatcagccatga