

FCGR3A-Fc fragment of IgG, low affinity IIIa, receptor (CD16a) Gene

Reference: PTXBC033678

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD16A; FCG3; FCGR3; FCGR3III; FCR-10; FCRIII; FCRIIIA; IGFR3; IMD20; low affinity immunoglobulin gamma Fc region receptor III-A; CD16a antigen; Fc fragment of IgG low affinity IIIa receptor; Fc fragment of IgG, low affinity III, receptor for (CD16); Fc fra		

For research use only.

GENE DESCRIPTION	Fc fragment of IgG, low affinity IIIa, receptor (CD16a)
GENE ID	2214
ORIGIN SPECIES	Human
ACCESSIONS	BC033678
VECTOR	pDONR223
NCBI SYMBOL	FCGR3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggcagctgctcctccaactgctctgctactttagttcagctggcatgcggactgaagatctccaaaggctgtggtgttctctggagcctcaatgggtacaggggtgctcgagaaggacagtgctgactctgaagtgccaggaggcctactccccctgaggacaattccacacagtggtttc
acaafgagagcctcatctcaagccaggcctcgagctacttattgacgctgccacagttgacgacagtgaggagtagagtgccagacaaacctctccacctcagtgacccgggtgcagctagaagtcctatcggctggctgtgtctccaggccccctcgggtggtgttcaaggagggaag
accctattcacctgaggtgtcacagctggaagaacactgctctgataaggtcacatatttacagaatgggcaaaggcagggaagtattttatcataattctgacttctacattccaaaggccacactcaagacagcggctcctacttctcagggggctgtgtgggagtaaaatgtgtcttcag
agactgtgaacatcacatcactcaaggtttgtcagtgcaacatctcatctttccacctgggtaccaagtcttctctgtgtgtggtgactccttttgcagtgacacagagactatatttctgtgaagacaacattcgaagctcaacaagagactggaaggaccataaatttaaatg
gagaaaggaccctcaagacaaatga