

C17orf66-chromosome 17 open reading frame 66 Gene

Reference: PTXBC033734

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf66; protein HEATR9; HEAT repeat-containing protein 9; HEAT repeat containing 9		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 66
GENE ID	256957
ORIGIN SPECIES	Human
ACCESSIONS	BC033734
VECTOR	pDONR223
NCBI SYMBOL	C17orf66
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctatgaaatcaactgatatctctgatgtctccaggtcaatgttctgtacccatggctggaatatccagacaagccaagaactcagaaaagccatggctcctgttctatctgcccttgcctgctaccagatgccaaaggaagagttcccccaagtcagagtgctggaggcagc
atccgagcaagccaactcagtcctgctacttcaagaaacctgagatctacacgcactggcagcagctgtatgatcagcagaggaagggagcctgagaagatgttggagaaatgagagatgactgttaggtacatcaaaaggatcatcaaacccatcaaaatgttccatct
cccaatgagcaagctgactataaaatctgagatgcgatccaggcccttagagcctacccaggacccccctgaatggcgaagattaagggaactcacaacaaagcctggaatctccagagagagatgagcagttctatgcagcacaggctctgggatgcttacgcatcagtgacaagtgtg
atggaggcactacagcaggtggtgagctgagaactgtggaggtcaggaaagaaggaacagggtctgagatgtgtggcatgtgggcacagacgaggtccccgggtgctag