

CHD9-chromodomain helicase DNA binding protein 9 Gene

Reference: PTXBC033770

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent helicase CHD9; AD013; CHD-9; CReMM; KISH2; PRIC320; chromodomain-helicase-DNA-binding protein 9; PPAR-alpha-interacting complex protein 320 kDa; PPAR{gamma}-interacting cofactor 320 kDa; chromatin remodeling factor CHROM1; chromatin-related		

For research use only.

GENE DESCRIPTION	chromodomain helicase DNA binding protein 9
GENE ID	80205
ORIGIN SPECIES	Human
ACCESSIONS	BC033770
VECTOR	pDONR223
NCBI SYMBOL	CHD9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttaataatcttttgtagccagctgaacatgtgttatctccacactctcagtttaattgtttccaatccatccccaaaacccaacccaatggtttgtttccagatgtatcagatggcagtcgaatgtggggccatcagacagctactaccatttcaaatcaaaatggatctctttccaccaaa
ggacattcacactctatgcataaaataaaagctttgtggcacaccatgactttgacctttatcaggccaatgaacaacaacacagtgacttctactacgctcacacaaaacagaaataatctcaaccaggcgagaattctcttagccagtcataaaattttatga