



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF38-ring finger protein 38 Gene

Reference: PTXBC033786

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

E3 ubiquitin-protein ligase RNF38; ring finger protein 38

For research use only.

GENE DESCRIPTION	ring finger protein 38
GENE ID	152006
ORIGIN SPECIES	Human
ACCESSIONS	BC033786
VECTOR	pDONR223
NCBI SYMBOL	RNF38
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr