



Reference: PTXBC033799

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCF-2; HCF2; host cell factor 2; C2 factor; host cell factor C2		
For research use only.			

GENE DESCRIPTION	host cell factor C2
GENE ID	29915
ORIGIN SPECIES	Human
ACCESSIONS	BC033799
VECTOR	pDONR223
NCBI SYMBOL	HCFC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr