



Reference: PTXBC033816

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATRC3; CAT-3; CAT3; cationic amino acid transporter 3; solute carrier family 7 (cationic amino acid transporter, y+ system), member 3; solute carrier family 7 member 3		

For research use only.

GENE DESCRIPTION	solute carrier family 7 (cationic amino acid transporter, y ⁺ system), member 3
GENE ID	84889
ORIGIN SPECIES	Human
ACCESSIONS	BC033816
VECTOR	pDONR223
NCBI SYMBOL	SLC7A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr