



Reference: PTXBC033821

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf44; kelch-like protein 36; kelch-like 36; kelch like family member 36		
For research use only.			

GENE DESCRIPTION	kelch-like 36 (Drosophila)
GENE ID	79786
ORIGIN SPECIES	Human
ACCESSIONS	BC033821
VECTOR	pDONR223
NCBI SYMBOL	KLHL36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]