

NIPBL-Nipped-B homolog (Drosophila) Gene

Reference: PTXBC033847

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIPBL, cohesin loading factor; CDLS; CDLS1; IDN3; IDN3-B; Scc2; nipped-B-like protein; Nipped-B homolog; SCC2 homolog; delangin; sister chromatid cohesion 2 homolog		

For research use only.

GENE DESCRIPTION	Nipped-B homolog (Drosophila)
GENE ID	25836
ORIGIN SPECIES	Human
ACCESSIONS	BC033847
VECTOR	pDONR223
NCBI SYMBOL	NIPBL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatgttgcagaaaatcagctccttaatcgaattgcaaatgtgtccagggtattttattactctcatgttaaaacacatttgaagaatcttggatttctgatagtaaaatcagaagtactctcatctgaatctgcaaaagtatatgataaagcgataaaccgaaaacaggagtt
cattttcatccaaacaaacactggacttctgaggagtgacatggttaattccaaatcacagaagaggtgaaaaggagtatatgaaaacagfatctagatttcaacttctcatggaacatctggaccctgatgaagaagaagaagaggagggtttcagctagcacaaatgctcggaac
aaagcaattacctcactgcttggaggaggcagccctaaaaataacagcagcagagacagaagatgatgaaagtgtggggaggatagaggaggaggcacttcaggggtgaggcggaggaggagtcacgtatttcgagcgtattacgtaa