

TCP10L-t-complex 10 (mouse)-like Gene

Reference: PTXBC033855

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf77; PRED77; TCP10A-2; T-complex protein 10A homolog 2; T-complex 10A-2; T-complex protein 10A-2; TCP10-like; t-complex 10 (a murine tcp homolog)-like; t-complex 10-like		

For research use only.

GENE DESCRIPTION	t-complex 10 (mouse)-like
GENE ID	140290
ORIGIN SPECIES	Human
ACCESSIONS	BC033855
VECTOR	pDONR223
NCBI SYMBOL	TCP10L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggcaggtcaactcaggccagggacccaaggaggcaccaccagaggaccgtgccaggagctggggctgtcatggagaagacagctgtggcagccgaggtctcacggaggactgcaacactggggagatgccaccattacagcagatcatcagactccac
caagagcttgggagacagaagtctctgtggctgatgttcatggaaaactccggagtcataatagctcttggaggagcagaacatggagctccgagaaaagctgagagctctgcagctgcagcgggtggaaagccaggaaagaatctgcagctccccacacgcggggcagaatcg
cacactctggcattggaacctgctttggaaaaattcacctctgtcagctgatgaagagacaatacccaatacgtggccacaagaatcagagtggcactctcctgggacaaagatcgtcatcaacaattcagctcctccaagccaatgagtttaagatagaagaattagctgtgga
aaacaccaccacaggaataagataaaaaatttccaggagacgtcaagacagaagacaacactactgggaaggccaactccctgtgcagagagacggggggtgtctga