

FCER1G-Fc fragment of IgE, high affinity I, receptor for; gamma polypeptide Gene

Reference: PTXBC033872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FCRG; high affinity immunoglobulin epsilon receptor subunit gamma; Fc epsilon receptor Ig; Fc fragment of IgE, high affinity I, receptor for; gamma polypeptide; Fc receptor gamma-chain; fc-epsilon RI-gamma; fcRgamma; fceRI gamma; immunoglobulin E receptor		

For research use only.

GENE DESCRIPTION	Fc fragment of IgE, high affinity I, receptor for; gamma polypeptide
GENE ID	2207
ORIGIN SPECIES	Human
ACCESSIONS	BC033872
VECTOR	pDONR223
NCBI SYMBOL	FCER1G
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattccagcagtggtctgtcttactccttttggtgaacaagcagcgccctggagagcctcagctctgtatatcctggtgcatcctgtttctgtatggaattgtcctcaccctctctactgtcgaactgaaggtatccaagtgcgaaagcagctataaccagctatgagaatcag
atggtgtttacacggcctgagcaccaggaaccaggagacttacgagactctgaagcatgagaaaccaccacagtag