

BST2-bone marrow stromal cell antigen 2 Gene

Reference: PTXBC033873

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD317; TETHERIN; bone marrow stromal antigen 2; BST-2; HM1.24 antigen; NPC-A-7; bone marrow stromal cell antigen 2		

For research use only.

GENE DESCRIPTION	bone marrow stromal cell antigen 2
GENE ID	684
ORIGIN SPECIES	Human
ACCESSIONS	BC033873
VECTOR	pDONR223
NCBI SYMBOL	BST2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatctactctgtatgactattgcagagtgcccatggaagacggggataagcgtgttaagcttctgctggggataggaattctgtgtctctgatcatcgtgattctgggggtgcccttgattatctcaccatcaaggccaacagcagggcctgccgggacggccttcgggcagtgtg
gagtgctgcaatgtcaccatctctgcaacaagagctgaccgagggccagaaagggcttcaggatgtggaggcccgccacctgcaaccacactgtgatggccctaattggcttccttgatgcagagaaggcccaaggacaaaagaaagtggaggagcttgaggagagat
cactacataaaccataagcttcaggacgctctgcagaggtggagcgactgagaagagaaaacaggctcttaagcgtgagaatcgcggacaagaagtactacccagctcccaggactccagctccgctgcggcgccccagctgctgattgtgctgctgggcctcagcgtctgctgc
agtga