

SPeM1-spermatid maturation 1 Gene

Reference: PTXBC033882

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf83; spermatid maturation protein 1; spermatid maturation 1		
For research use only.			

GENE DESCRIPTION	spermatid maturation 1
GENE ID	374768
ORIGIN SPECIES	Human
ACCESSIONS	BC033882
VECTOR	pDONR223
NCBI SYMBOL	SPeM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccatggttgagcggccgagggccgagtgaggcctcgtatcacaactgcaacagcaacagctgccaggacctgggcaactctgtctgttgctgtgggcctcatctgcattaacattagcatcaatagtgacctgctctggagccgattccgtggtgtcttataccaagtgttcc
atgataccatttggagaagaagctcccaagtcattactcagaaagcagaccagcccccctaaagaagcagagttctctcgtcagtcctctcgggtgcaccatggacctgtgatgatgctgtgtcccccagctcaccgccatcgcctcagagctctccacacgctgtgctc
actgccagtagcttggcctcctgacactgatgacgagaagcctcatcagtaaccagcctatctgtctcctaccactgggatgtccctgaggactgggaagcctccaacacactcaggggacctgggttccctggtctcagatgccccagagtcctcccccagaccatccgttccagcc
taccgttagaggaaagccctcaaacagcctatggtcggagctgggcccctaaagggcctatgtgtatctgtgaacccccacctccagccctgaggctcctagccacaagaacgggtggggaggggggcggtgccagaggcagaggcggtcagtaaccagcctgtccagctccca
ccctggccagcagtcacccctgaatttcccgccaccgctcctcagggccgaatagtgatgtgcccgggacatgagacggcggttcgggaactgacctgggaggtggaggccctgtccggtgtacccctagcctctggtatccagactgccgaggagacaagaagaattg
ggtgtaccgttccttaactgggaggtga