



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DACT3-dapper, antagonist of beta-catenin, homolog 3 (*Xenopus laevis*) Gene

Reference: PTXBC034052

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

DAPPER3; RRR1; dapper homolog 3; antagonist of beta-catenin Dapper homolog 3; arginine-rich region 1 protein; dapper, antagonist of beta-catenin, homolog 3; dishevelled binding antagonist of beta catenin 3

For research use only.

GENE DESCRIPTION	dapper, antagonist of beta-catenin, homolog 3 (<i>Xenopus laevis</i>)
GENE ID	147906
ORIGIN SPECIES	Human
ACCESSIONS	BC034052
VECTOR	pDONR223
NCBI SYMBOL	DACT3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccggggctctctgctccgggtgagccctgagccggggccgctgctggggctgctgagggtagtacctggccgggctctgcgagtacattggctccgggagaggcaggagtagccgcgtgcagcaggccgtcggctgcccagcccggaaatggggggccggaggccga
ggcagaggaggagcccgatgaagatgaagatcgcggcgccggccgcccagcggccctcgaggagcagctggagccctgctggtctctctgggacctggggacagcagctgggagacctgagcctggagagctggggcctggaacaggagagaggcggtgtagctc
ggcgctctatgaagatccagctctacaggaggtccagattcaccacccctcaacctctgctgggacagctgggcttctggatccagctctatgctgcctcggtgctctgagccccggcatctatgcaggtgagagggcccaagtcctatgtaaagggtgagggtggggccctg
ggcgccggggcaagggcagctg

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1