

TOMM5-translocase of outer mitochondrial membrane 5 homolog (yeast) Gene

Reference: PTXBC034247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf105; Tom5; bA613M10.3; mitochondrial import receptor subunit TOM5 homolog; mitochondrial outer membrane protein TOM5; translocase of outer mitochondrial membrane 5 homolog; translocase of outer mitochondrial membrane 5		

For research use only.

GENE DESCRIPTION	translocase of outer mitochondrial membrane 5 homolog (yeast)
GENE ID	401505
ORIGIN SPECIES	Human
ACCESSIONS	BC034247
VECTOR	pDONR223
NCBI SYMBOL	TOMM5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcggattgagggcctcgccgaagctggacccggaggagatgaaacgggaagatgcgcgaggatgtgatctctccatacggaaactttctcatctacgtggccctcctgcgagtcactccattatcttaagaaattggacagcatatga