

FMR1NB-fragile X mental retardation 1 neighbor Gene

Reference: PTXBC034320

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT37; NY-SAR-35; NYSAR35; fragile X mental retardation 1 neighbor protein; cancer/testis antigen 37; sarcoma antigen NY-SAR-35; fragile X mental retardation 1 neighbor		

For research use only.

GENE DESCRIPTION	fragile X mental retardation 1 neighbor
GENE ID	158521
ORIGIN SPECIES	Human
ACCESSIONS	BC034320
VECTOR	pDONR223
NCBI SYMBOL	FMR1NB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcacataggagaaagcgaaggagggaataggagaagtcaccgtgccatgcgtgtggtcacttagagctggcaacttatgagttggcggcaactgagtcgaatcccagagcagccatcctgatacagagccgcatggtgacaggcctcagccaggatggcgggaa
tctctaaagatgcgggtcagcaaacctttgggatgctcatgctccatttggatcctgctgttcgtgtgtactactgtctcctgctcgggtcctcatattttgtgcttgcgaatggacatactgcccacagtgaaaatgctcatgccaatctctggaagaagattccgcatg
gaagctttgctgaatttttttccacaacttgaatcaggggaaatcaggtggcaagccttgaatgagctgcaagatcttagtgagagtgaaatgtttgagacacaaatgctgttttcatcatcggggaccacagctcaaatgttttctccattagagatgtgcctaaacagatgatg
caaatgtttggccttggtgcgatcagccttatcctggtatgtgtcccatatttgcgcgtctcttttctggaggagcgaaccggccgatgattacaaggcaggacacagagttgtaacgggtttgaagaacaaagaagggaagcgaagagggaagtctgaaatgttacagaaagcagca
agaggacgtgaggaaacatggtgacgtag