



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ELOVL3-elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 3
Gene

Reference: PTXBC034344

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3-keto acyl-CoA synthase ELOVL3; CIG-30; CIG30; elongation of very long chain fatty acids protein 3; ELOVL FA elongase 3; cold-inducible glycoprotein of 30 kDa; elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 3; very long chai		

For research use only.

GENE DESCRIPTION	elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 3
GENE ID	83401
ORIGIN SPECIES	Human
ACCESSIONS	BC034344
VECTOR	pDONR223
NCBI SYMBOL	ELOVL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcacagccatgaatgtctcacatgaagtaaatcagctgttcacgccccttaactctgagctgtccaaagacatgagcccttttcgaggagatattgggcaacctctatccccatagccctgatctacctgtgtctcatcgtctgtgggcagaaactcatgaaggaacgcagggtctcaa
cctgcaagggcctctatcctctgtgctcttctcctgtcaatctcagatctctgtgggcagtgagagtggtgggcattatggggacgtgtgctacttaccggggccctaaagcaaacctgtgtgcttcaacaactcatcgataattccacagtcacaaattctgtctcgtgtcttctctcagcaag
gtcatagaactcgcgagacacagcttcatcatctctgctgaagcgcgcacatcttattactctgtaccacacagcagcagctgtctgtctacaaagcttttgatatacaagaacaaagctgcctgcaggaagcgtgttctgtcatcatgaactctgtgtgttcatgccatcatgacacctctaacct
ctgaaggtgcgcacagtggaagcccccaagatctgcccctgatctcatcaccagcctgcagatctgtgttgatagctcgtacagctcctcagctacgtcagctgacagagatgccacaccacgatggaacattattctgtgtcttcatctgttatatgaccttattcat
ctgtctgtccacttctctgcacacacatcagtcgcccactgtcgaagcgaagacagagagccagctga

More info: Please contact us for any question: contact@proteogenix.fr