

ARL14-ADP-ribosylation factor-like 14 Gene

Reference: PTXBC034354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARF7; ADP-ribosylation factor-like protein 14; ADP-ribosylation factor 7; ADP-ribosylation factor-like 14; ADP ribosylation factor like GTPase 14		

For research use only.

GENE DESCRIPTION	ADP-ribosylation factor-like 14
GENE ID	80117
ORIGIN SPECIES	Human
ACCESSIONS	BC034354
VECTOR	pDONR223
NCBI SYMBOL	ARL14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggttcgctgggttctaaatccgcaaaccaagccaagttctcttttgggacttgactcagctgggaagtctactctctttataaattaaagcttgctaaagatattaccaccatccctacaatagggttcaatgtgaaatgatcaggttgaaaggaattttcactcacagtctgg
gatgttgaggacagagaaaaatgagaactgtttgggctgttactgtgagaacaccgatgggctgtgtatgtgtggacagtacagacaaacagcagactggaagagttctcagagacagtttgagcacatttgaagaatgaacacattaaatgtcctgtgttctattagccaacaaac
aagacatgcctggagctctgactgtgagacatcaccagaatgttcaaatgaagaagcttgcagtgaccggaactggtatgtgcaaccctgtgtgcccacagggggaggggtggccaggggttcaggaaattaactggattgtgaagagccacatgaatcaagaggagac
actttggcgttctcaagcagaactga